EQUIPMENT CALIBRATION, STANDARDIZATION OR CHECK RECORD

Revised 12/96, 3/11, 3/14

				PROCEDURE #40
Date:		Checke	d by:	
Equipment: MECHANICAL SHAKER (Thoroughness) T 27				
Previous check date:		Next due:	Freque	ncy: 12 Months Max
Mechanical Shaker serial number:				
Scale readable to 0.1 g, 20,000 g capacity #:				
Action recommended:	Repair	Replace No	one	
MINUTES SHAKEN				
Sieve Serial No.	Sieve Size	(B) Wt. Pass hand shake	(A) Total Wt.	0.5% Maximum (B / A)*100